

 Dentallegierung	Dichte	Wachsumrechnungstabelle Wachsgewicht x Dichte + 10 %																															
	g/cm <sup>3</sup>	Wachsgewicht in Gramm																															
		0,4	0,5	0,6	0,7	0,8	0,9	1	1,1	1,2	1,3	1,4	1,5	1,6	1,7	1,8	1,9	2	2,2	2,4	2,6	2,8	3	3,2	3,4	3,6	3,8	4	4,2	4,4	4,6	4,8	5
<b>Goldguss-Legierungen für Kronen, Brücken und Goldfüllungen</b>																																	
Wegold Biopur	19,3	8	11	13	15	17	19	21	23	25	28	30	32	34	36	38	40	42	47	51	55	59	64	68	72	76	81	85	89	93	98	102	106
Wegold i	17,7	8	10	12	14	16	18	19	21	23	25	27	29	31	33	35	37	39	43	47	51	55	58	62	66	70	74	78	82	86	90	93	97
Wegold JN	15,4	7	8	10	12	14	15	17	19	20	22	24	25	27	29	30	32	34	37	41	44	47	51	54	58	61	64	68	71	75	78	81	85
Wegold Biocomp	16,7	7	9	11	13	15	17	18	20	22	24	26	28	29	31	33	35	37	40	44	48	51	55	59	62	66	70	73	77	81	85	88	92
Wegold M-PF	15,6	7	9	10	12	14	15	17	19	21	22	24	26	27	29	31	33	34	38	41	45	48	51	55	58	62	65	69	72	76	79	82	86
Wegold B-SG	15,4	7	8	10	12	14	15	17	19	20	22	24	25	27	29	30	32	34	37	41	44	47	51	54	58	61	64	68	71	75	78	81	85
Wegold M	15,8	7	9	10	12	14	16	17	19	21	23	24	26	28	30	31	33	35	38	42	45	49	52	56	59	63	66	70	73	76	80	83	87
Wegold G	13,8	6	8	9	11	12	14	15	17	18	20	21	23	24	26	27	29	30	33	36	39	43	46	49	52	55	58	61	64	67	70	73	76
<b>Universal-Legierungen für Spezialkeramik (z.B. cosmica)</b>																																	
Wegold Bio-Activity	16,9	7	9	11	13	15	17	19	20	22	24	26	28	30	32	33	35	37	41	45	48	52	56	59	63	67	71	74	78	82	86	89	93
Wegold Copynorm	16,7	7	9	11	13	15	17	18	20	22	24	26	28	29	31	33	35	37	40	44	48	51	55	59	62	66	70	73	77	81	85	88	92
Wegold Norm Supra	15,8	7	9	10	12	14	16	17	19	21	23	24	26	28	30	31	33	35	38	42	45	49	52	56	59	63	66	70	73	76	80	83	87
Wegold Bio-EVO	16,6	7	9	11	13	15	16	18	20	22	24	26	27	29	31	33	35	37	40	44	47	51	55	58	62	66	69	73	77	80	84	88	91
Wegold Bioconcept	16,6	7	9	11	13	15	16	18	20	22	24	26	27	29	31	33	35	37	40	44	47	51	55	58	62	66	69	73	77	80	84	88	91
Wegold Norm	16,1	7	9	11	12	14	16	18	19	21	23	25	27	28	30	32	34	35	39	43	46	50	53	57	60	64	67	71	74	78	81	85	89
Wegold M-NF	15,6	7	9	10	12	14	15	17	19	21	22	24	26	27	29	31	33	34	38	41	45	48	51	55	58	62	65	69	72	76	79	82	86
Wegold Biomedic PDF	15,6	7	9	10	12	14	15	17	19	21	22	24	26	27	29	31	33	34	38	41	45	48	51	55	58	62	65	69	72	76	79	82	86
Wegold NF IV	13,7	6	8	9	11	12	14	15	17	18	20	21	23	24	26	27	29	30	33	36	39	42	45	48	51	54	57	60	63	66	69	72	75
Wegold PM-NF	10,6	5	6	7	8	9	10	12	13	14	15	16	17	19	20	21	22	23	26	28	30	33	35	37	40	42	44	47	49	51	54	56	58
<b>Gold-Platin-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
Wegold GS	19,1	8	11	13	15	17	19	21	23	25	27	29	32	34	36	38	40	42	46	50	55	59	63	67	71	76	80	84	88	92	97	101	105
Wegold Bio-Attraction	19,0	8	10	13	15	17	19	21	23	25	27	29	31	33	36	38	40	42	46	50	54	59	63	67	71	75	79	84	88	92	96	100	105
Wegold Bio-REVO	18,9	8	10	12	15	17	19	21	23	25	27	29	31	33	35	37	40	42	46	50	54	58	62	67	71	75	79	83	87	91	96	100	104
Wegold Bio-Compose	18,9	8	10	12	15	17	19	21	23	25	27	29	31	33	35	37	40	42	46	50	54	58	62	67	71	75	79	83	87	91	96	100	104
Wegold H	18,2	8	10	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40	44	48	52	56	60	64	68	72	76	80	84	88	92	96	100
Wegold U	17,7	8	10	12	14	16	18	19	21	23	25	27	29	31	33	35	37	39	43	47	51	55	58	62	66	70	74	78	82	86	90	93	97
Wegold HG plus	18,6	8	10	12	14	16	18	20	23	25	27	29	31	33	35	37	39	41	45	49	53	57	61	65	70	74	78	82	86	90	94	98	102
<b>Palladium-Basis-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
Wegold N2	12,1	5	7	8	9	11	12	13	15	16	17	19	20	21	23	24	25	27	29	32	35	37	40	43	45	48	51	53	56	59	61	64	67
Wegold DG	11,6	5	6	8	9	10	11	13	14	15	17	18	19	20	22	23	24	26	28	31	33	36	38	41	43	46	48	51	54	56	59	61	64
Wegold MT	11,5	5	6	8	9	10	11	13	14	15	16	18	19	20	22	23	24	25	28	30	33	35	38	40	43	46	48	51	53	56	58	61	63
Wegold Ag	11,3	5	6	7	9	10	11	12	14	15	16	17	19	20	21	22	24	25	27	30	32	35	37	40	42	45	47	50	52	55	57	60	62
<b>Ecosystem-Legierungen für Spezialkeramik (z.B. cosmica)</b>																																	
Wegold econorm	13,9	6	8	9	11	12	14	15	17	18	20	21	23	24	26	28	29	31	34	37	40	43	46	49	52	55	58	61	64	67	70	73	76
Wegold EcoBeta	12,8	6	7	8	10	11	13	14	15	17	18	20	21	23	24	25	27	28	31	34	37	39	42	45	48	51	54	56	59	62	65	68	70
Wegold ecogold	12,7	6	7	8	10	11	13	14	15	17	18	20	21	22	24	25	27	28	31	34	36	39	42	45	47	50	53	56	59	61	64	67	70
Wegold Ecologic	12,7	6	7	8	10	11	13	14	15	17	18	20	21	22	24	25	27	28	31	34	36	39	42	45	47	50	53	56	59	61	64	67	70
Wegold ecoflex	12,6	6	7	8	10	11	12	14	15	17	18	19	21	22	24	25	26	28	30	33	36	39	42	44	47	50	53	55	58	61	64	67	69
Wegold Ecotec	10,7	5	6	7	8	9	11	12	13	14	15	16	18	19	20	21	22	24	26	28	31	33	35	38	40	42	45	47	49	52	54	56	59
Wegold ecoplus	10,6	5	6	7	8	9	10	12	13	14	15	16	17	19	20	21	22	23	26	28	30	33	35	37	40	42	44	47	49	51	54	56	58
Wegold ecolight	10,7	5	6	7	8	9	11	12	13	14	15	16	18	19	20	21	22	24	26	28	31	33	35	38	40	42	45	47	49	52	54	56	59
<b>Ecosystem-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
Wegold eco-connect	12,1	5	7	8	9	11	12	13	15	16	17	19	20	21	23	24	25	27	29	32	35	37	40	43	45	48	51	53	56	59	61	64	67
Wegold ecosave	11,3	5	6	7	9	10	11	12	14	15	16	17	19	20	21	22	24	25	27	30	32	35	37	40	42	45	47	50	52	55	57	60	62

		0,4	0,5	0,6	0,7	0,8	0,9	1	1,1	1,2	1,3	1,4	1,5	1,6	1,7	1,8	1,9	2	2,2	2,4	2,6	2,8	3	3,2	3,4	3,6	3,8	4	4,2	4,4	4,6	4,8	5
<b>Goldguss-Legierungen für Kronen, Brücken und Goldfüllungen</b>																																	
Aurofluid® 2 PF	15,9	7	9	10	12	14	16	17	19	21	23	24	26	28	30	31	33	35	38	42	45	49	52	56	59	63	66	70	73	77	80	84	87
Pontor® MPF	15,5	7	9	10	12	14	15	17	19	20	22	24	26	27	29	31	32	34	38	41	44	48	51	55	58	61	65	68	72	75	78	82	85
Neocast® 3	15,5	7	9	10	12	14	15	17	19	20	22	24	26	27	29	31	32	34	38	41	44	48	51	55	58	61	65	68	72	75	78	82	85
Protor® 3	15,0	7	8	10	12	13	15	17	18	20	21	23	25	26	28	30	31	33	36	40	43	46	50	53	56	59	63	66	69	73	76	79	83
Pontor® 2	14,2	6	8	9	11	12	14	16	17	19	20	22	23	25	27	28	30	31	34	37	41	44	47	50	53	56	59	62	66	69	72	75	78
Yellow Special	12,5	6	7	8	10	11	12	14	15	17	18	19	21	22	23	25	26	28	30	33	36	39	41	44	47	50	52	55	58	61	63	66	69
<b>Gold-Platin-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
BioEthic®	18,9	8	10	12	15	17	19	21	23	25	27	29	31	33	35	37	40	42	46	50	54	58	62	67	71	75	79	83	87	91	96	100	104
V-Gnathos® Supra	19,1	8	11	13	15	17	19	21	23	25	27	29	32	34	36	38	40	42	46	50	55	59	63	67	71	76	80	84	88	92	97	101	105
V-Gnathos® Plus	18,8	8	10	12	14	17	19	21	23	25	27	29	31	33	35	37	39	41	45	50	54	58	62	66	70	74	79	83	87	91	95	99	103
Esteticor® Royal H	18,6	8	10	12	14	16	18	20	23	25	27	29	31	33	35	37	39	41	45	49	53	57	61	65	70	74	78	82	86	90	94	98	102
Esteticor® Avenir®	18,7	8	10	12	14	16	19	21	23	25	27	29	31	33	35	37	39	41	45	49	53	58	62	66	70	74	78	82	86	91	95	99	103
V 92	18,0	8	10	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40	44	48	51	55	59	63	67	71	75	79	83	87	91	95	99
Esteticor® Cosmor H	18,2	8	10	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40	44	48	52	56	60	64	68	72	76	80	84	88	92	96	100
<b>Gold-Palladium-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
V-Classic	16,4	7	9	11	13	14	16	18	20	22	23	25	27	29	31	32	34	36	40	43	47	51	54	58	61	65	69	72	76	79	83	87	90
V-Delta® SF	14,5	6	8	10	11	13	14	16	18	19	21	22	24	26	27	29	30	32	35	38	41	45	48	51	54	57	61	64	67	70	73	77	80
Esteticor® Accurate 40	13,6	6	7	9	10	12	13	15	16	18	19	21	22	24	25	27	28	30	33	36	39	42	45	48	51	54	57	60	63	66	69	72	75
<b>Palladium-Basis-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
Cerapall® 6	11,7	5	6	8	9	10	12	13	14	15	17	18	19	21	22	23	24	26	28	31	33	36	39	41	44	46	49	51	54	57	59	62	64
Ceradelta® 2	11,3	5	6	7	9	10	11	12	14	15	16	17	19	20	21	22	24	25	27	30	32	35	37	40	42	45	47	50	52	55	57	60	62
<b>Implant-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
Esteticor® Implant® 76	16,9	7	9	11	13	15	17	19	20	22	24	26	28	30	32	33	35	37	41	45	48	52	56	59	63	67	71	74	78	82	86	89	93
Esteticor® Implant® 58	15,1	7	8	10	12	13	15	17	18	20	22	23	25	27	28	30	32	33	37	40	43	47	50	53	56	60	63	66	70	73	76	80	83
Esteticor® Implant® 32	13,1	6	7	9	10	12	13	14	16	17	19	20	22	23	24	26	27	29	32	35	37	40	43	46	49	52	55	58	61	63	66	69	72
<b>Ecosystem-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
Esteticor® Economic	14,2	6	8	9	11	12	14	16	17	19	20	22	23	25	27	28	30	31	34	37	41	44	47	50	53	56	59	62	66	69	72	75	78
Esteticor® Blancor	11,3	5	6	7	9	10	11	12	14	15	16	17	19	20	21	22	24	25	27	30	32	35	37	40	42	45	47	50	52	55	57	60	62
EcoDelta	11,3	5	6	7	9	10	11	12	14	15	16	17	19	20	21	22	24	25	27	30	32	35	37	40	42	45	47	50	52	55	57	60	62
Esteticor® Actual	11,2	5	6	7	9	10	11	12	14	15	16	17	18	20	21	22	23	25	27	30	32	34	37	39	42	44	47	49	52	54	57	59	62



Dentallegierung

Dichte

g/cm<sup>3</sup>

Wachsumrechnungstabelle

Wachsgewicht x Dichte + 10 %

Wachsgewicht in Gramm

		0,4	0,5	0,6	0,7	0,8	0,9	1	1,1	1,2	1,3	1,4	1,5	1,6	1,7	1,8	1,9	2	2,2	2,4	2,6	2,8	3	3,2	3,4	3,6	3,8	4	4,2	4,4	4,6	4,8	5
<b>Goldguss-Legierungen für Kronen, Brücken und Goldfüllungen</b>																																	
Multiceram PDF	16,7	7	9	11	13	15	17	18	20	22	24	26	28	29	31	33	35	37	40	44	48	51	55	59	62	66	70	73	77	81	85	88	92
Aurotec SG	15,4	7	8	10	12	14	15	17	19	20	22	24	25	27	29	30	32	34	37	41	44	47	51	54	58	61	64	68	71	75	78	81	85
Aurotec M	15,8	7	9	10	12	14	16	17	19	21	23	24	26	28	30	31	33	35	38	42	45	49	52	56	59	63	66	70	73	76	80	83	87
<b>Universal-Legierungen für Spezialkeramik (z.B. cosmica)</b>																																	
Multiceram Norm	16,7	7	9	11	13	15	17	18	20	22	24	26	28	29	31	33	35	37	40	44	48	51	55	59	62	66	70	73	77	81	85	88	92
Multiceram Supra	15,8	7	9	10	12	14	16	17	19	21	23	24	26	28	30	31	33	35	38	42	45	49	52	56	59	63	66	70	73	76	80	83	87
Multiceram Best	13,9	6	8	9	11	12	14	15	17	18	20	21	23	24	26	28	29	31	34	37	40	43	46	49	52	55	58	61	64	67	70	73	76
Aurotec NF IV	13,7	6	8	9	11	12	14	15	17	18	20	21	23	24	26	27	29	30	33	36	39	42	45	48	51	54	57	60	63	66	69	72	75
Ecoceram Gold	12,2	5	7	8	9	11	12	13	15	16	17	19	20	21	23	24	25	27	30	32	35	38	40	43	46	48	51	54	56	59	62	64	67
Ecoceram Ag	10,6	5	6	7	8	9	10	12	13	14	15	16	17	19	20	21	22	23	26	28	30	33	35	37	40	42	44	47	49	51	54	56	58
<b>Gold-Platin-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
Bioceratec N	19,0	8	10	13	15	17	19	21	23	25	27	29	31	33	36	38	40	42	46	50	54	59	63	67	71	75	79	84	88	92	96	100	105
Ceratec H	18,2	8	10	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40	44	48	52	56	60	64	68	72	76	80	84	88	92	96	100
Ceratec U	17,7	8	10	12	14	16	18	19	21	23	25	27	29	31	33	35	37	39	43	47	51	55	58	62	66	70	74	78	82	86	90	93	97
<b>Palladium-Basis-Legierungen für Metallkeramik (z.B. classica, unica)</b>																																	
Ecoceram DG	11,6	5	6	8	9	10	11	13	14	15	17	18	19	20	22	23	24	26	28	31	33	36	38	41	43	46	48	51	54	56	59	61	64

Verarbeitungsanleitungen  
finden Sie auf [www.wegold.de](http://www.wegold.de)

Wegold Edelmetalle GmbH

Nibelungenstr. 5  
90530 Wendelstein

Telefon +49 (0)129 4030-0  
Telefax +49 (0)129 4030-40

info@wegold.de • www.wegold.de

QMf 4.5-8 Rev. f vom 11.01.2019